


MUN.PZL




Basic kit with n.6 tips
for osteotomy, osteoplasty, extraction



MUN.PZ - Classic version
(no led)



MUN.PZL - Led version



SPECIFICATIONS

- ▶ Newest generation piezo handpiece with high luminosity LED light
- ▶ Electronic control with accurate adjustment of functional parameters
- ▶ Extended service life of parts that require autoclave sterilization
- ▶ Practical peristaltic pump with fast tube replacement and high efficiency
- ▶ Easy touch keyboard and wide display

Operating frequency piezo: 27 kHz + 37 kHz automatic scan.
Piezo handpiece power: 60 W.
Main functions piezo mode: n. 3 programs, control of power and irrigation, led activation.

CARATTERISTICHE

- ▶ Manipolo piezo di ultima generazione con luce LED ad alta intensità
- ▶ Controllo elettronico con precisa regolazione dei parametri funzionali
- ▶ Sterilizzazione: lunga durata delle parti autoclavabili
- ▶ Pratica pompa peristaltica con cambio tubo rapido e massima efficienza
- ▶ Tastiera touch e ampio display

Frequenza di funzionamento piezo: 27 kHz + 37 kHz scansione automatica.
Potenza utile al manipolo piezo: 60 W.
Funzioni principali modalità piezo: n. 3 programmi, regolazione potenza e irrigazione, gestione led.

The High Quality Made in Italy

MARIOTTI & C. SRL

Implant-surgery motors and piezoelectric - manufacturer since 1973
Via Seganti 73 - 47121 Forlì - Italy
Tel +39 0543 781811
info@mariotti-italy.com
www.mariotti-italy.com

